



Associated Course Details

Name: Toy Cleaning

ID: 070403012

Area of Focus / Skillset: Toys

Supplies

- Acetone
- Scent free dish detergent
- Dry towels
- Enzymatic cleaner
- Small, soft microfiber wash rag or disposable rag
- Sealable plastic bag (larger than toy itself)

Tips / Notes / Warnings

The best way to remove odors is to prevent the odors from infiltrating toys. A general light clean-out and an immediate wash after soiling helps to reduce odor absorption.

NOTE: Regular use of acetone may result in the loss of the shiny luster of a toy.

WARNING: Do NOT use acetone or bleach on PVC, latex, rubber, low grade silicone, and ABS plastics.

For additional information, please visit:
<https://www.squarepegtoys.com/support/how-do-i-remove-the-odor-from-my-plug/>

Steps / Procedures

To remove odors from silicone toys:

1. Wash toy thoroughly using soap and water to remove any lube and debris from the play session. Dry completely.
2. Wet a soft, disposable rag with acetone and gently rub affected areas.
3. Rinse the toy and reassess for odors.
4. If the odor persists, treat with enzymatic cleaner.
 - For small toys, fully immerse toy in small container of enzymatic solution for 30 minutes.
 - For large toys,
 - Place toy inside sealable plastic bag.
 - Spray affected areas with an enzymatic cleaner.
 - Seal the bag and wait four hours.
 - Reassess for odors and repeat if necessary.
5. Wash a final time with soap and water and rinse thoroughly to remove any residual cleaners.
6. Remove from wet areas and dry thoroughly before storing.